

# **BL21 Star<sup>™</sup> (DE3) One Shot<sup>®</sup>** **BL21 Star<sup>™</sup> (DE3)pLysS One Shot<sup>®</sup>** **Chemically Competent Cells**

Catalog nos. C6010-03, C6020-03

Version D  
September 17, 2002  
25-0402



[www.invitrogen.com](http://www.invitrogen.com)  
[tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)



# Table of Contents

Table of Contents .....	iii
Kit Contents and Storage .....	v
<b>Introduction .....</b>	<b>1</b>
Overview.....	1
Choosing a BL21 Star™ Strain.....	3
<b>Methods .....</b>	<b>5</b>
Basic Transformation Procedure .....	5
Expression Guidelines.....	7
<b>Appendix .....</b>	<b>10</b>
Testing Transformation Efficiency .....	10
Technical Service.....	11
Purchaser Notification .....	13
Product Qualification .....	15
References .....	16



## Kit Contents and Storage

### Available Kits

The table below lists the One Shot® chemically competent cell kits covered by this manual. The transformation efficiency is calculated as number of transformants per µg of pUC19 plasmid DNA.

Item	Reactions	Transformation Efficiency	Catalog no.
BL21 Star™(DE3)	20	1 x 10 <sup>8</sup> cfu/µg	C6010-03
BL21 Star™(DE3)pLysS	20	1 x 10 <sup>8</sup> cfu/µg	C6020-03

### Shipping/Storage

Each One Shot® kit is shipped on dry ice. Upon receipt, store the kit at -80°C.

### Kit Contents

The table below describes the items included in each of the One Shot® chemically competent *E. coli* kits described above. **Store at -80°C.**

Item	Composition	Amount
SOC Medium (may be stored at room temperature or +4°C)	2% tryptone 0.5% yeast extract 10 mM NaCl 2.5 mM KCl 10 mM MgCl <sub>2</sub> 10 mM MgSO <sub>4</sub> 20 mM glucose	6 ml
Chemically competent cells	--	21 x 50 µl
pUC19 Control DNA	10 pg/µl in 5 mM Tris-HCl, 0.5 mM EDTA, pH 8	50 µl

### Genotypes

**BL21 Star™(DE3):** F<sup>-</sup> *ompT hsdS<sub>B</sub> (r<sub>B</sub><sup>-</sup>m<sub>B</sub><sup>-</sup>) gal dcm rne131* (DE3)

**BL21 Star™(DE3)pLysS:** F<sup>-</sup> *ompT hsdS<sub>B</sub> (r<sub>B</sub><sup>-</sup>m<sub>B</sub><sup>-</sup>) gal dcm rne131* (DE3) pLysS (Cam<sup>R</sup>)

The DE3 designation means the strains contain the λ DE3 lysogen which carries the gene for T7 RNA polymerase under control of the *lacUV5* promoter. IPTG is required to induce expression of the T7 RNA polymerase.

The two strains carry a mutated *rne* gene (*rne131*) which encodes a truncated RNase E enzyme that lacks the ability to degrade mRNA, resulting in an increase in mRNA stability (see page 1).

The two strains are *E. coli* B/r strains and do not contain the *lon* protease. They are also deficient in the outer membrane protease, OmpT. The lack of these proteases reduces degradation of heterologous proteins expressed in the strains.

BL21 Star™(DE3)pLysS carries the pLysS plasmid which produces T7 lysozyme (see page 1). The BL21 Star™(DE3) strain does **not** carry a plasmid expressing T7 lysozyme.



# Introduction

## Overview

---

### Introduction

The BL21 Star™(DE3) and BL21 Star™(DE3)pLysS *E. coli* strains are derived from the BL21(DE3) and BL21(DE3)pLysS strains, respectively. In addition to the λ DE3 lysogen which allows high-level expression of T7-regulated genes, the BL21 Star™ strains also contain the *rne131* mutation to enhance the expression capabilities of each strain.

---

### *rne131*

The *rne* gene encodes the RNase E enzyme, an essential, 1061 amino acid *E. coli* endonuclease which is involved in rRNA maturation and mRNA degradation as a component of a protein complex known as a “degradosome” (Grunberg-Manago, 1999; Lopez *et al.*, 1999). Various studies have shown that the N-terminal portion of RNase E (approximately 584 amino acids) is required for rRNA processing and cell growth while the C-terminal portion of the enzyme (approximately 477 amino acids) is required for mRNA degradation (Kido *et al.*, 1996; Lopez *et al.*, 1999). The *rne131* mutation (present in the BL21 Star™ strains) encodes a truncated RNase E which lacks the C-terminal 477 amino acids of the enzyme required for mRNA degradation (Kido *et al.*, 1996; Lopez *et al.*, 1999). Thus, mRNAs expressed in the RNase E-defective BL21 Star™ strains exhibit increased stability. When heterologous genes are expressed in these strains from T7-based expression vectors, the yields of recombinant proteins generally increase.

---

### pLysS

The pLysS plasmid carried by the BL21 Star™(DE3)pLysS strain produces T7 lysozyme to reduce basal level expression of the gene of interest. pLysS confers resistance to chloramphenicol (Cam<sup>R</sup>) and contains the p15A origin. This origin allows pLysS to be compatible with pUC- or pBR322-derived plasmids.

---

### Expressing Heterologous Genes

**The BL21 Star™(DE3) and BL21 Star™(DE3)pLysS strains are suitable for high-level recombinant protein expression.** Due to the increase in stability of mRNAs, we have observed higher basal expression of heterologous genes in BL21 Star™ strains than in some BL21 strains; therefore, these strains may not be useful for expression of toxic genes. To choose an appropriate BL21 Star™ strain to use for expression of your gene of interest, see pages 2 and 3.

**Note:** Basal expression levels of heterologous genes are generally higher in BL21 Star™(DE3) cells than in BL21 Star™(DE3)pLysS cells. However, the overall yield of recombinant protein is generally higher in BL21 Star™(DE3) cells than in BL21 Star™(DE3)pLysS cells.

---

*continued on next page*

## Overview, continued

---



The BL21 Star™ strains are useful to express heterologous genes from any T7-based expression vector. In general, we recommend using the BL21 Star™(DE3) strain to express heterologous genes from low-copy number, T7-based plasmids (*e.g.* pET vectors). We have observed toxicity when expressing some heterologous genes from high-copy number plasmids (*e.g.* Invitrogen's pCR®T7 vectors) in BL21 Star™(DE3) cells. These effects are alleviated when BL21 Star™(DE3)pLysS cells are used; therefore, we recommend using the BL21 Star™(DE3)pLysS strain to express heterologous genes from high-copy number, T7-based plasmids.

---

### General Handling

Be extremely gentle when working with competent cells. Competent cells are highly sensitive to changes in temperature or mechanical lysis caused by pipetting. Transformation should be started immediately after thawing the cells on ice, and all mixing should be done by swirling or tapping the tube gently, not by pipetting.

---



### Important


BL21 Star™(DE3) and BL21 Star™(DE3)pLysS cells require IPTG to induce expression of the T7 RNA polymerase from the *lacUV5* promoter. IPTG is available from Invitrogen (Catalog no. 15529-019). For more information, see our Web site or call Technical Service (see page 11).

---

# Choosing a BL21 Star™ Strain

## Choosing a BL21 Star™ Strain for Expression

The flowchart below provides some general guidelines to help you choose the appropriate BL21 Star™ strain to use in your expression experiments.



For more information about indications of toxicity in bacterial cells, see the **Expression Guidelines**, page 8.

*continued on next page*

## Choosing a BL21 Star™ Strain, continued

---

### Comparison with Other BL21 Strains

We have compared expression of many heterologous genes from T7-based expression vectors in the RNase E-defective BL21 Star™ strains (BL21 Star™(DE3) or BL21 Star™(DE3)pLysS) to non-RNase E-defective BL21 strains (BL21(DE3) or BL21(DE3)pLysS). In most cases, the yield of recombinant protein obtained from the BL21 Star™ strain ranged from 2 to 10-fold greater than that from the corresponding non-RNase E-defective BL21 strain. The table below provides a representative list of recombinant proteins which we have found to be expressed to higher levels in the BL21 Star™ strains when compared to the non-RNase E-defective BL21 strains. Note that the pET-based vectors contain the *lacI* gene and the T7/*lac* promoter.

Recombinant Protein	Expression Vector	BL21 Star™ Strain	BL21 Strain
Firefly luciferase	pET-based	BL21 Star™(DE3)	BL21(DE3)
Vaccinia topoisomerase I	pET-based	BL21 Star™(DE3)	BL21(DE3)
<i>E. coli</i> β-galactosidase ( <i>lacZ</i> )	pCR®T7-based	BL21 Star™(DE3)pLysS	BL21(DE3)pLysS
Human RARγ	pCR®T7-based	BL21 Star™(DE3)pLysS	BL21(DE3)pLysS
Human TFIIB ( <i>GTF2B</i> )	pCR®T7-based	BL21 Star™(DE3)pLysS	BL21(DE3)pLysS
Human BDNF receptor ( <i>NTRK2</i> )	pCR®T7-based	BL21 Star™(DE3)pLysS	BL21(DE3)pLysS
Human Cdk5 ( <i>CDK5</i> )	pCR®T7-based	BL21 Star™(DE3)pLysS	BL21(DE3)pLysS

**Note:** Some heterologous genes may not be expressed to higher levels in the BL21 Star™ strains when compared to non-RNase E-defective BL21 strains.

---

# Methods

## Basic Transformation Procedure

---

### Introduction

A basic transformation protocol for BL21 Star™(DE3) and BL21 Star™(DE3)pLysS cells is provided below. Once you have selected transformants, we recommend that you proceed directly to expression using your own protocol. **Note that BL21 Star™(DE3) and BL21 Star™(DE3)pLysS are designed to be used for expression, not cloning or subcloning.**

---

### Materials Supplied by the User

- Plasmid DNA (ready for transformation)
  - 42°C water bath
  - 37°C shaking and non-shaking incubator
  - Ice bucket with ice
  - Spectrophotometer or similar device to measure optical density of the bacterial cultures
  - Microcentrifuge tube rack (optional)
- 

### Before Starting

- Prepare LB agar plates containing the appropriate concentration of antibiotic (to select for your plasmid). If you are transforming DNA into BL21 Star™(DE3)pLysS cells, add 34 µg/ml chloramphenicol (to select for pLysS)
  - Equilibrate a water bath to 42°C
  - Warm the vial of SOC medium to room temperature
  - Place the plates in a 37°C incubator to remove excess moisture (use two plates for each transformation)
- 

*continued on next page*

## Basic Transformation Procedure, continued

---

### Basic Transformation Procedure

1. Thaw on ice one vial of One Shot<sup>®</sup> cells per transformation.
  2. Add 5 to 10 ng of DNA, in a volume of 1 to 5  $\mu$ l, to the cells and mix by tapping gently. **Do not mix cells by pipetting.**
  3. Incubate the vial(s) on ice for 30 minutes.
  4. Incubate for exactly 30 seconds in the 42°C water bath. **Do not mix or shake.**
  5. Remove vial(s) from the 42°C bath and quickly place on ice.
  6. Add **250  $\mu$ l** of pre-warmed SOC medium to the vial(s). **Note:** SOC is a rich medium; practice good sterile technique to avoid contamination.
  7. Place the vial(s) in a microcentrifuge rack and secure with tape. Place the rack on its side in a shaking incubator and shake the vial(s) at 37°C for 1 hour at 225 rpm.
  8. Plate 20 to 200  $\mu$ l each of the transformation reaction onto two LB plates containing the appropriate antibiotic and 34  $\mu$ g/ml chloramphenicol (if using BL21 Star<sup>™</sup> (DE3) pLysS). Plate two different volumes to ensure well-spaced colonies on at least one plate. The remaining transformation reaction may be stored at +4°C and plated out the next day, if needed.
  9. Invert the plates and incubate at 37°C overnight.
  10. Select a transformant and culture as described on page 8. **Note:** Expression can vary between clones. You may wish to characterize additional transformants.
-

# Expression Guidelines

---

## Introduction

If you have an expression protocol for the plasmid that you are working with, we recommend that you use your own protocol. This section provides some general guidelines for the use of T7 RNA polymerase-based expression plasmids in BL21 Star™(DE3) or BL21 Star™(DE3)pLysS cells.

---



Transform your expression plasmid into a strain that does not bear the gene for T7 RNA polymerase (*i.e.* TOP10, DH5α) and maintain your construct in this strain. Use BL21 Star™(DE3) or BL21 Star™(DE3)pLysS cells **for expression only**.

---

## BL21 Star™ Strains

The BL21 Star™ strains are suitable for high-level recombinant protein expression. In general, the yield of recombinant protein obtained from BL21 Star™(DE3) cells is higher than the yield from BL21 Star™(DE3)pLysS cells.

---

## T7 RNA Polymerase and Toxic Genes

In the BL21 Star™(DE3) and BL21 Star™(DE3)pLysS strains, the T7 polymerase gene is controlled by the IPTG-inducible *lacUV5* promoter. Because of the extremely high activity of T7 RNA polymerase and the increased stability of mRNAs, some basal level expression of the gene of interest will likely occur in uninduced cells. This creates problems in cases where the gene of interest is toxic to bacterial cells. In these cases, expression of the toxic gene under uninduced conditions leads to selection of cells that express the lowest levels of the toxic gene. These cells are often unable to express high levels of the gene of interest upon IPTG induction of the T7 polymerase.

The BL21 Star™(DE3)pLysS strain produces T7 lysozyme which helps to reduce basal levels of T7 RNA polymerase. Although levels are reduced, the cells may still contain a small amount of T7 RNA polymerase.

---

## Using BL21 Star™(DE3)pLysS

We recommend using the BL21 Star™(DE3)pLysS strain if:

- You are using a high-copy, T7-based vector to express your gene of interest
  - You observe growth inhibitory effects (*i.e.* toxicity) when using BL21 Star™(DE3) (see the next page)
  - You are expressing a known toxic gene
- 

*continued on next page*

## Expression Guidelines, continued

---

### Expression Guidelines

The following guidelines assume that expression of your gene is not toxic to *E. coli*. If you are working with a toxic gene, some additional precautions may be taken (see the next page).

1. Following transformation (see page 6), pick 3 or 4 transformants and culture them in 5 ml LB medium containing the appropriate antibiotic to select for your expression plasmid. If you are using BL21 Star™(DE3)pLysS, remember to add 34 µg/ml chloramphenicol to select for pLysS. Grow at 37°C with shaking until the OD<sub>600</sub> reaches 0.6 to 1.0.
2. Use these cultures to inoculate fresh LB medium containing the appropriate antibiotic to an OD<sub>600</sub> of 0.05-0.1 (~1:20 dilution of the initial culture). This dilution allows the cells to quickly return to logarithmic growth and reach the appropriate cell density. Use a volume appropriate for taking time points, if desired.

**Note:** If you are using BL21 Star™(DE3)pLysS, you do not need to include chloramphenicol in these cultures. Generally, the cells will not lose the pLysS plasmid during the limited number of cell doublings that occur in the growth and induction stages.

3. Use the remainder of the initial culture to create a glycerol stock of the clones. Once you have identified the clone that best expresses your protein, you can use the glycerol stock to perform additional expression experiments.
4. Grow the cultures until they reach mid-log (OD<sub>600</sub>~0.4, 2 to 3 hours).
5. Induce the cultures by adding IPTG to a final concentration of 0.5 mM and culture for an additional 2-3 hours. You may also take time points to analyze for optimal expression of your protein.
6. Analyze clones by western blot or enzymatic assay to determine which clone best expresses your protein of interest. Use the glycerol stock created from this clone for expression experiments (see above). If you find that expression levels in subsequent inductions decrease, or you find that you lose your plasmid, your protein may be toxic to *E. coli* (see below for additional information).

---

### Indications of Toxicity

When expressing recombinant proteins in BL21 Star™ strains, we generally assume that the recombinant protein is toxic to bacterial cells when any of the following occurs:

- No transformants are obtained after following Steps 1-9 of the **Basic Transformation Protocol** on page 6 **OR** a combination of large and small, irregular colonies appears on the plate
- The initial culture does not grow (see Step 1 of the **Expression Guidelines**, above)
- It takes longer than 5 hours after a 1:20 dilution of the initial culture for the fresh culture to reach an OD<sub>600</sub>=0.4 (see Steps 2 and 4 of the **Expression Guidelines**, above)
- The cells lyse after induction with IPTG (see Step 5 of the **Expression Guidelines**, above)

---

*continued on next page*

## Expression Guidelines, continued

---

### Precautions

Listed below are several precautions that may be taken to prevent problems resulting from basal level expression of a toxic gene of interest. The guidelines below assume that you are using BL21 Star™(DE3)pLysS for expression. These precautions are designed to ensure that BL21 Star™(DE3)pLysS cells bearing the gene of interest are cultured for as short a time as possible before IPTG induction. These methods all assume that the T7 expression plasmid has been correctly designed and created.

- Propagate and maintain your expression plasmid in a strain that does not contain T7 RNA polymerase (*i.e.* TOP10, DH5 $\alpha$ , etc.).
  - Perform a fresh transformation of BL21 Star™(DE3)pLysS cells before each induction experiment.
  - Following transformation of BL21 Star™(DE3)pLysS cells using the **Basic Transformation Protocol** on page 6, pick 3 or 4 transformants directly into fresh LB medium containing the appropriate antibiotic. When the culture reaches OD<sub>600</sub>=0.4, induce expression of the recombinant protein by adding IPTG to a final concentration of 0.5 mM.
- 

### Transformation/ Expression Protocol for Toxic Genes

This alternative protocol may be used with BL21 Star™(DE3)pLysS cells. Other protocols are suitable, depending on your needs.

#### Transformation

1. Follow the basic transformation protocol on page 6 through Step 7.
2. After outgrowth of the transformation reaction (page 6, Step 7), add the **entire** transformation reaction (300  $\mu$ l) to 50 to 200 ml of LB containing 100  $\mu$ g/ml ampicillin and 34  $\mu$ g/ml chloramphenicol, prewarmed to 37°C.

#### Induction

3. Incubate the cells with shaking at 37°C until the cells reach mid-log phase (OD<sub>600</sub>=0.3). **Note:** Doubling times may vary (30 to 90 minutes) depending on the protein expressed.
  4. Add IPTG to a final concentration of 0.5 to 1 mM and grow for 2-3 more hours. You may take time points, if desired.
  5. Harvest cells by centrifugation and use immediately or store the cell pellet at -80°C.
- 

### Other Alternatives

If you are using BL21 Star™(DE3)pLysS cells and observe significant toxicity, you may want to try using the BL21-AI™ strain (Catalog no. C6070-03) available from Invitrogen to express your recombinant protein of interest. The BL21-AI™ strain contains a chromosomal insertion of the gene encoding T7 RNA polymerase into the *araB* locus of the *araBAD* operon, allowing expression of T7 RNA polymerase to be tightly regulated by L-arabinose. For more information about the BL21-AI™ strain, see our Web site ([www.invitrogen.com](http://www.invitrogen.com)) or call Technical Service (see page 11).

---

# Appendix

## Testing Transformation Efficiency

---

### Introduction

To test the transformation efficiency of the competent cells contained in the One Shot<sup>®</sup> kit, use the supercoiled pUC19 plasmid supplied with the kit as described below. An extra vial of cells is included for this purpose.

---

### Before Starting

- Prepare LB agar plates containing 50 µg/ml ampicillin.
  - Equilibrate a water bath to 42°C.
  - Warm the vial of SOC medium to room temperature.
  - Place the plates in a 37°C incubator to remove excess moisture (use two plates for each transformation).
- 

### Transformation

Follow the transformation protocol on page 6 to transform pUC19 into BL21 Star<sup>™</sup>(DE3) or BL21 Star<sup>™</sup>(DE3)pLysS. Use the specific modifications below.

- Transform cells with 1 µl (10 pg) of pUC19
  - Plate 50 µl each onto two LB plates containing 50 µg/ml ampicillin.
  - Calculate the transformation efficiency as transformants per µg of plasmid (see below). The cells should have an efficiency of  $1 \times 10^8$  transformants/µg of supercoiled plasmid.
- 

### Calculation

Use the formula below to calculate transformation efficiency.

$$\frac{\# \text{ of colonies}}{10 \text{ pg transformed DNA}} \times \frac{10^6 \text{ pg}}{\mu\text{g}} \times \frac{300 \mu\text{l transformed cells}}{X \mu\text{l plated}} = \frac{\# \text{ transformants}}{\mu\text{g plasmid DNA}}$$

---

# Technical Service

---

## World Wide Web



Visit the [Invitrogen Web Resource](#) using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

**<http://www.invitrogen.com>**

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

---

## Contact Us

For more information or technical assistance, please call, write, fax, or email. Additional international offices are listed on our web page ([www.invitrogen.com](http://www.invitrogen.com)).

### Corporate Headquarters:

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008  
USA

Tel: 1 760 603 7200

Tel (Toll Free): 1 800 955 6288

Fax: 1 760 602 6500

E-mail: [tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)

### Japanese Headquarters:

Invitrogen Japan K.K.  
Nihonbashi Hama-Cho Park  
Bldg. 4F  
2-35-4, Hama-Cho, Nihonbashi

Tel: 81 3 3663 7972

Fax: 81 3 3663 8242

E-mail: [jpinfo@invitrogen.com](mailto:jpinfo@invitrogen.com)

### European Headquarters:

Invitrogen Ltd  
Inchinnan Business Park  
3 Fountain Drive  
Paisley PA4 9RF, UK

Tel: +44 (0) 141 814 6100

Tech Fax: +44 (0) 141 814 6117

E-mail: [eurotech@invitrogen.com](mailto:eurotech@invitrogen.com)

---

## MSDS Requests

To request an MSDS, visit our web site ([www.invitrogen.com](http://www.invitrogen.com)) and follow the instructions below.

1. On the home page, go to the left-hand column under 'Technical Resources' and select 'MSDS Requests'.
  2. Follow instructions on the page and fill out all the required fields.
  3. To request additional MSDSs, click the 'Add Another' button.
  4. All requests will be faxed unless another method is selected.
  5. When you are finished entering information, click the 'Submit' button. Your MSDS will be sent within 24 hours.
- 

*continued on next page*

## Technical Service, continued

---

### Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives.

Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

**Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

---

# Purchaser Notification

---

## Introduction

Use of any BL21 Star™ *E. coli* strain is covered under a number of different licenses including those detailed below.

---

## Information for European Customers

The BL21 Star™(DE3) and BL21 Star™(DE3)pLysS strains are genetically modified and carry the bacteriophage λ DE3 lysogen containing the T7 RNA polymerase gene. As a condition of sale, this product must be in accordance with all applicable local legislation and guidelines including EC Directive 90/219/EEC on the contained use of genetically modified organisms.

---

## Limited Use Label License No. 30: T7 Expression System

The composition and/or use of this product is claimed in one or more patents U.S. Patent Nos. 4,952,496 and 5,693,489 and 5,869,320 licensed to Invitrogen Corporation, by Brookhaven Science Associates, LLC.

The T7 expression system is based on technology developed at Brookhaven National Laboratory under contract with the U.S. Department of Energy, and is the subject of patents and patent applications assigned to Brookhaven Science Associates, LLC (BSA,). By provisions of the Distribution License Agreement granted to Invitrogen covering said patents and patent applications, Invitrogen grants you a non-exclusive sub-license under patents assigned to BSA for the use of this technology, including the enclosed materials, based upon the following conditions: 1 – these materials are to be used for non-commercial research purposes only. A separate license under patents owned by BSA is required for any commercial use, including the use of these materials for research purposes or production purposes by any commercial entity. Information about commercial license may be obtained from The Office of Technology Transfer, Brookhaven National Laboratory, Bldg. 475D, P.O. Box 5000, Upton, New York 11973-5000. Phone (516) 344-7134. 2 – No materials that contain the cloned copy of the T7 gene 1, the gene for T7 RNA polymerase, may be distributed further to third parties outside of your laboratory, unless the recipient receives a copy of this sub-license and agrees to be bound by its terms. This limitation applies to strains BL21(DE3), BL21(DE3)pLysS and BL21(DE3)pLysE, CE6, BL21-SI Competent Cells and any derivatives that are made of them.

You may refuse this sub-license by returning this product unused in which case Invitrogen accept return of the products with a full refund. By keeping or using this product, you agree to be bound by the terms of this license.

---

*continued on next page*

## Purchaser Notification, continued

---

**Limited Use Label  
License No. 50:  
BL21 Star™**

Invitrogen Corporation (“Invitrogen”) has an exclusive license to sell the *rne131* genotype to scientists for **research or commercial evaluation purposes only**, under the terms described below. Use of the *rne131* genotype by commercial entities for Commercial Purposes, as defined below, beyond evaluation requires the user to obtain a commercial license as detailed below. Before using the BL21 Star™ strain, please read the following license agreement. If you do not agree to be bound by its terms, contact Invitrogen within 10 days for authorization to return the unused BL21 Star™ product and to receive a full credit.

The *rne131* genotype is covered by one or more French patents or patent applications and corresponding foreign patents or patent applications owned by CNRS. Information about commercial licenses may be obtained from: Centre National de la Recherche Scientifique (CNRS), 3, rue Michel Ange, 75794 PARIS Cedex 16, France, REF: L00084, Contact: [Dale.Roche@fist.fr](mailto:Dale.Roche@fist.fr)

Invitrogen grants you a non-exclusive license to use the enclosed BL21 Star™ strain for research or for commercial evaluation purposes only. The strain is being transferred to you in furtherance of, and reliance on, such license. You may not use the strain, or the materials contained therein, for any Commercial Purpose, as defined below, beyond a one-year evaluation without a license for such purpose from CNRS.

If you are a commercial entity, each of your laboratories is allowed a one-year evaluation (not free use) period, after which time this right automatically terminates. To use any portion of BL21 Star™ strain for a Commercial Purpose as defined below, commercial entities must obtain a commercial license from CNRS for each of their laboratories. Contact information for commercial entities purchasing a BL21 Star™ strain will be provided to CNRS, who may contact them during the evaluation period regarding their desire for a commercial license. Commercial Purposes include: Any use of *rne131* in a Commercial Product; Any use of *rne131* in the manufacture of a Commercial Product; Any sale of *rne131* or products having used in a commercial process the *rne131* genotype.

Access to the BL21 Star™ strain must be limited solely to those officers, employees and students of your institution who need access to perform the above-described research or evaluation. You must inform each such officer, employee and student of the provisions of this license agreement and require them to agree to be bound by the provisions of this license agreement. You may not distribute the BL21 Star™ strain to others, even to those within your own institution. You may only transfer modified, altered, or original material from the cell line to a third party following written notification of, and written approval from, Invitrogen so that the recipient can be licensed. You may not assign, sub-license, rent, lease or otherwise transfer this license agreement or any of the rights or obligation there under, except as expressly permitted by Invitrogen and CNRS.

This license agreement is effective until terminated. You may terminate it at any time by destroying all BL21 Star™ products in your control. It will also terminate automatically if you fail to comply with the terms and conditions of the license agreement. You shall, upon termination of the license agreement, destroy all BL21 Star™ strains in your control, and so notify Invitrogen in writing.

---

# Product Qualification

---

## Product Specifications

The following criteria are used to qualify One Shot® BL21 Star™(DE3) and BL21 Star™(DE3)pLysS cells:

50 µl of competent cells are transformed with 10 pg of supercoiled pUC19 plasmid DNA. Transformed cultures are plated on LB plates containing 50 µg/ml ampicillin and the transformation efficiency is calculated. Test transformations are performed in triplicate. Transformation efficiency should be:

- Greater than  $1 \times 10^8$  cfu/µg DNA for BL21 Star™(DE3) cells
- Greater than  $1 \times 10^8$  cfu/µg DNA for BL21 Star™(DE3)pLysS cells

Untransformed cells are plated on:

- LB plates containing 50 µg/ml ampicillin to verify the absence of ampicillin resistant contamination.
  - LB plates as a lawn to verify the absence of phage contamination.
  - LB plates containing 34 µg/ml chloramphenicol for selection of pLysS (for BL21 Star™(DE3)pLysS only)
-

## References

---

Grunberg-Manago, M. (1999). Messenger RNA Stability and its Role in Control of Gene Expression in Bacteria and Phages. *Annu. Rev. Genet.* 33, 193-227.

Kido, M., Yamanaka, K., Mitani, T., Niki, H., Ogura, T., and Hiraga, S. (1996). RNase E Polypeptides Lacking a Carboxyl-terminal Half Suppress a *mukB* mutation in *Escherichia coli*. *J. Bacteriol.* 178, 3917-3925.

Lopez, P. J., Marchand, I., Joyce, S. A., and Dreyfus, M. (1999). The C-terminal Half of RNase E, Which Organizes the *Escherichia coli* Degradosome, Participates in mRNA Degradation but not rRNA Processing *in vivo*. *Mol. Microbiol.* 33, 188-199.

---

©2000-2002 Invitrogen Corporation. All rights reserved.